# TAURUS Transport stretcher

1

10

PROMA REHA®



## TAURUS

The TAURUS transport trolley is designed for fast and safe transport of bariatric patients weighing up to 285 kg.

The ergonomic and intuitive controls allow immediate use of all functions. The handrails are equipped with the SAFE GAPS safety system, which eliminates the risk of pinching or trapping the hands during folding. The large castors and lightweight design ensure high mobility and easy handling.



#### High load capacity 320 kg

The solid and robust design with a load capacity of up to 320 kg creates a secure support for patient care. Very high durability is ensured thanks to the materials and technologies used.

### 5th wheel

The fifth wheel ensures easy and intuitive handling of the stretcher. This feature is particularly useful in rugged areas and for handling by one person.

## Optimal width for bariatric patients

The width of the stretcher guarantees comfortable transport even for larger patients and at the same time allows perfect access for the medical staff to the patient.

#### Comfortable positioning

Back positioning is aided by a gas piston rod that allows instant adjustment to the desired position without physical effort.



#### Safe side rails

The Side rails along the patient surface can be easily folded down to provide immediate access to the patient. Thanks to the SAFE GAPS safety system, the risk of pinching the arm during folding is eliminated.



#### Universal accessory holders



Fixation belt set

#### Infusion stand holder



#### Oxygen bottle holder



#### Easy height adjustment

The mattress platform lift is controlled by an easy-to-use foot pedal accessible from both sides.





Mattress from medical leatherette The comfortable mattress is upholstered in medical leatherette, which has anti-static protection with antimicrobial properties.

Folding transport handle at the head The transport handling handle behind the patient's head allows for quick folding when needed.

#### TENTE brand castors

High-quality anti-static castors with centering brake, also suitable for outdoor use, reduce the physical effort of the staff and, thanks to the covers, facilitate cleaning and significantly extend their service life.

#### Plastic chassis cover

The plastic cover makes daily cleaning of the stretcher easier and faster. It also protects the chassis mechanisms, thus extending the life of the individual components of stretcher.



#### Paper roll holder

The hygiene paper roll holder can be found under the mattress platform.

Bumper wheels The bumper wheels located at the corners of the mattress platform are a great help especially during transport.









#### Technical parameters

Outer dimensions	
Mattress platform	
Mattress platform height	57 - 94 cm
Lift of mattress platform	Hydraulic
Back part positioning	0 - 85°
Safe operating load	
Maximum patient weight	285 kg
Weight	110 kg
Type of side rails	Foldable
Castors	Metal, 200 mm
Type of brake	Central

Material of mattress platform	Medical leatherette
Set of restraint straps	Yes (3 x)
5th wheel	Yes
Bumper wheels	Yes
IV pole holders	Yes (2 x)
Trasport holder of IV pole under the mattress platform	Yes
Transport handrail foldable to head	Yes
Transport handrail to foot	Yes
Oxygen bottle holder	Yes
Horizontal accessories holder	Yes
Paper roll holder	Yes

PROMA REHA®



www.promareha.com

Riegrova 342, 552 03 Česká Skalice, Czech Republic Tel: +420 775 776 122, info@promareha.com