



URGENT

Transport stretcher for fast and safe transport of patients

Medical staff will appreciate the easy handling, ergonomic and intuitive controls and the full range of functions. Large castors and a lightweight design ensure high mobility and easy handling. The handrails are equipped with the SAFE GAPS safety system, which eliminates the risk of pinching or trapping the hands during folding. URGENT can be complemented with a wide range of optional accessories.



High load capacity

The solid and robust construction with a load capacity of up to 230 kg forms a secure support for patient care.





Infusion stand transport holder below the mattress platform

The infusion stand holder is conveniently placed below the mattress platform so that it is always immediately available.

Choice of colour

The choice of colour design for the leatherette of the couch allows patients and staff to create an aesthetically beautiful and inviting space that can help improve mood, comfort and reduce stress.







Hydraulic height adjustment

For easier access to the patient, the couch is equipped with a hydraulic lift of the patient area ergonomically located in the middle of the couch.

Safe side rails

The protective side rails along the mattress platform can be easily folded down to provide immediate access to the patient. Thanks to the SAFE GAPS safety system, the risk of pinching or trapping the arm during folding is eliminated.





Medical leatherette mattress platform The comfortable mattress is upholstered in medical leatherette, which has anti-static protection with antimicrobial properties.

> Foldable head transport handle The transport handle behind the patient's head allows quick folding when necessary, which is especially appreciated during resuscitation.

TENTE brand castors

High-quality medical castors with central brake, reducing physical effort for staff. Castors are made of permanent non-marking material. Castor covers facilitate cleaning and extend the life of the castors.

> Oxygen bottle holder The oxygen bottle can be safely stored and transported thanks to the special holder below the mattress platform.

Plastic chassis cover

The plastic cover makes daily cleaning easier and faster. It also protects the chassis mechanisms, thus extending the life components and the overall reliability of the stretcher.

Bumper wheels

The bumper wheels located at the corners of the mattress platform are a great help especially during transport.













Technical parameters

Outer dimensions	
Dimensions of mattress platform	
Mattress platform positioning	Hydraulic
Backrest positioning	
Backrest mechanism	Piston rod
Safe operating load	230 kg
Maximum patient weight	195 kg
Weight	90 kg
Types of side rails	Folding, hinged, stainless steel
Castors	Metal, 200 mm
Brake type	Central

Mattress platform material	Upholstered medical leatherette
Fixation belt set	Yes
Transport handle at head	Yes foldable
Horizontal accessories holder	Yes
Transport brackets for infusion stand under the mattress platform	
Oxygen bottle holder on chassis	Yes
Mattress with anti-static cover	Yes
Bumper wheels	Yes
Plastic chassis cover	Optional
Transport handle at foot	Optional
Hygiene paper roll holder	Optional

PROMA REHA®



www.promareha.com

Riegrova 342, 552 03 Česká Skalice, Czech Republic Tel: +420 775 776 122, info@promareha.com