



SALVATOR


Bed for standard nursing care.

SALVATOR electric care beds are specially designed for standard nurceing care, especially for people with reduced mobility or those who need assistance with everyday activities. These beds are equipped with various features such as electric positioning, robust folding side rails and a controller with colour-coded functions that make the nurses job easier and help the clients to be independent.



Increased load capacity

The solid and stable construction with a safe operating load of up to 250 kg provides safe support even when caring for patients with higher weight. Thanks to the materials and technologies used, extremely high durability is ensured even under repeated increased loads.



Extension of the mattress platform


Extending the mattress platform increases comfort even for taller patients and makes it easier for staff to handle the patient.


Width of the mattress platform

The optimal width of the mattress platform ensures maximum comfort for the patient and easy handling of the bed.









Comfortable positioning

The increased range of the bed lift reduces the strain on staff and facilitates the client handling. The reduced minimum height of the mattress platform assists the client in exiting the bed. The electric positioning of the back and thigh section in combination with the mechanical positioning of the calf section allows a wide range of adjustment of the individual parts of the mattress platform. The staff and client can conveniently adjust the individual positions using the hand control.




REHA


Foldable side rails

The safety of the client on the bed is ensured by folding side rails. The easy folding of the side rails with two independent movements significantly reduces the risk of accidental lowering. The folded side rails are positioned sufficiently below the mattress level to allow easy exit from the bed.

Universal accessories holders



Bumper wheels







Plastic undercarriage cover

Plastic wheels with central brake

Bed linen shelf



Choice of design and decor

Being able to choose the design and colour decor of care beds allows you to create an esthetic and welcoming space that can help improve mood, comfort and reduce stress.





Technical parameters

Mattress platform (length x width)	
Outer bed dimensions (length x width)	
Mattress platform height	
Mattress platform extension	20 cm (optional)
Weight	118 kg
Safe working load	250 kg
Maximum patient weight	215 kg
CPR of the back	Ano
The bed is designed for application environ	ments
Bedding aream	etal slats (plastic removable optionally)
Types of bedheads	UNI P, UNI K

Bed equipment

Castors	Metal, plastic (optional)
Bumper wheels	Yes (optional)
Central brake	Yes (optional)
Universal accessory holder	Yes 2x (optional)
Plastic hooks for universal accessory holder	Yes 6x (optional)
Chassis cover	Yes (optional)
Accessory holders	Yes 4x
Bed linen shelf	Yes (optional)
Patient control	Yes
Backup battery	Yes (optional)

PROMA REHA®

www.promareha.com

Riegrova 342, 552 03 Česká Skalice, Czech Republic Tel: +420 775 776 122, info@promareha.com

04/2025-a-EN