



DOMUS 3

width 88 - 90 cm

III. deg. pressure ulcers

load capacity 140 - 200 kg



height 10,2 - 20,3 cm

air mattress

Dynamic anti-decubitus air mattress



Domus 3 combines alternating and static therapy modes to improve pressure relief care and reduce caregiver workload. It includes a low-pressure alert for patient safety and a weight indicator dial for easy adjustment. Additionally, it effectively wicks away moisture from the surface to aid in microclimate management, leading to a more comfortable experience and better prevention of pressure injury without losing cell firmness. In 2015, a clinical case study¹ involving 10 patients demonstrated that the Domus 3 system effectively improved moisture lesions in affected patients.



TECHNICAL PARAMETERS

- active system up to stage III decubitus and for patients with a high risk of developing decubitus
- possibility of adjusting the patient's weight
- the system is placed on a foam mattress base, but contains a pocket for a foam insert, thanks to which it can be placed directly on the bed frame
- CPR valve for quick deflation of the mattress
- static mode, transport mode and low pressure alarm (light)
- 17 air cells with pressure change after 10 minutes, 3 cells under the head without alternation
- cover 4-way elastic, vapor-permeable, antibacterial with 360° zipper and washable at 90°C

MATTRESS

size 190 - 210 x 85 - 90 cm
height 10,2 - 20,3 cm
weight 3,1 - 6,9 kg
type overlay / replacement
load capacity 140 - 200 kg
top cover material Nylon / PU
cell material Nylon / TPU

COMPRESSOR

dimensions 29 x 18,5 x 12,6 cm
mode alternate / static
noise <35 dB
holder integrated
weight 2,2 kg
material Fire Retardant ABS

