





The BABY newborn cot is the ideal solution for ensuring the closest contact between mother and baby in the first days of life. The bed is equipped with a transparent cot tray that allows you to watch your baby and at the same time provides him with plenty of natural light, which is important for his healthy development.

Height adjustment

The height adjustment of the cot allows parents to easily adjust the height of the cot to suit their needs.





Removable transparent cot tray

Thanks to the transparent cot tray, not only will you be able to see the baby, but also the sun's rays, so necessary for every baby. The plastic cot tray is made of a health-safe material that fully complies with the highest requirements and standards for use in contact with the smallest babies. There are holes in the bottom of the cot tray to ensure optimal air circulation under the mattress. The cot tray can be easily removed by pulling it out of the cot structure.

Trendelenburg / Antitrendelenburg

Sleeping in the right position is very important for newborns in the first months of life. The inclination of the bath minimises the risk of reflux and keeps the airway clear. The baby can sleep peacefully, which contributes to his overall health and comfort.

±15°

Handle for easy handling For easy handling and safe moving, the cot

is equipped with a handling handle.

 PROMA REHA

Castors

The braked castors make the cot easy to transport and safe to handle. So you can always have your baby with you no matter what you do.



Choice of design

The possibility of choosing the colour design of the cot allows you to create an aesthetic and pleasant space that can help to improve the mood, comfort and reduce stress not only for the baby but also for his mother during the stay in the maternity ward.







Technical parameters

Suitable for newborns	up to 3 months of age
External dimensions (I x w x h)	109 x 51 x 95 - 125 cm
Mattress platform dimensions (I x w)	75 x 40 cm
Height of mattress platform	69 - 99 cm
Mattress dimensions	74 x 39 x 4 cm
Height adjustability	Gas piston
Bed surface tilt control	Gas piston
Trendelenburg / Antitrendelenburg	15° / 15°
Safe operating load	10 kg
Castors	65 mm, plastic, braked
Cot tray	Removable, transparent
Weight	16 kg

PROMA REHA®



www.promareha.com

Riegrova 342, 552 03 Česká Skalice, Czech Republic Tel: +420 775 776 122, info@promareha.com

06/2023-a-EN